

Wheelchair to Bed Transfer

- 1) Bring w/c to foot of bed and lock w/c first
- 2) Turn wheels inward
- 3) Remove footrests
- 4) Make your feet are in front of the patient. Patients feet are touching the ground.
- 5) Ensure patient has on non-skid slippers or non-skid socks on
- 6) Stand in front of patient
- 7) Counting with patient 1,2,3 cueing patient.
- 8) Using your leg muscles lift and pivot patient to bed
- 9) Have patient sit on side of bed
- 10) With one hand unlock w/c and move it back, while the other arm is supporting the patient
- 11) With one arm supporting patient shoulders and other under the knees on the count of 1,2
3 using your leg muscles and one swift motion assist patient back into bed.
- 12) Removes shoes /slippers